Instructions For Use

EN. Do Not Reuse
DE. Nicht wiederverwenden
FR. Ne le réutilisez pas
IT. Non riutilizzare
ET. Mitte taaskasutada
PL. Nie używaj ponownie
DA. Må ikke genbruges
FI. Älä käytä uudelleen
NL. Niet hergebruiken
CS. Nepoužívejte znovu
LV. Neizmantoj atkārtoti
LT. Nenaudokite pakartotinai
PT. Não reutilize
SV. Återanvänd inte
SK. Nepoužívajte znova
TR. Yeniden kullanmayın
ES. No reutilizer
EL. Μην επαναχρησιμοποιείτε
HU. Ne használd újra

En. Caution
DE. Vorsicht
FR. Attention
IT. Attenzione
ET. Ettevaatust
PL. Zwóć uwagę
DA. Advarsel
FI. Varoitus
NL. Let op
CS. Upozornění
LV. Pievērš uzmanību
LT. Įspėjimas
PT. Cuidado
SV. Försiktighet
SK. Dávaj pozor
TR. Dikkat etmek
ES. Precaución
EL. Προσοχή
HU. Vigyázat

En. Use By
DE. Verwendung bis
FR. Utiliser par :
IT. Data di scadenza
ET. Kasutamine
PL. Do użycia przez
DA. Anvendelse af
FI. Käytä mennessä
NL. Te gebruiken door
CS. Použítí do
LV. Lietot līdz
LT. Naudoti iki
PT. Usado por
SV. Användning av
SK. Ak čeezte použiť
TR. Tarafından kullanılabacak
ES. Usado por
EL. Χρήση από
HU. Hasznája

EN. Guide Catheter/Minimum Inner Diameter
DE. Führungskatheter/Minimaler Innendurchmesser
FR. Cathéter-guide/diamètre intérieur minimum
IT. Catetere guida/Diametro interno minimo
ET. Juhtkateeter/Minimaalne sisemine läbimõõt
PL. Cewnik prowadzący/Minimalna średnica wewnętrzna
DA. Vejledende kateter/Minimal indre diameter
FI. Ohjauskatetri/Minimi sisähalkaisija
NL. Geleidingskatheter/Minimale binnendiameter
CS. Vodic íkat dr/Minimální vnitřní průměr
LV. Vadošais katetra/Minimālais iekšējais diametrs
LT. Kreipiamasis kateteris/Minimalus vidinis skersmuo
PT. Cateter guia/Diâmetro interno m ímimo
SV. Vägledande kateter/Minimal inre diameter
SK. Vodiáci katéter/Minimálny vnútorný priemer
TR. Kılavuz kateter/Minimum iç çap
ES. Catéter guía/Diámetro interior m ímimo
EL. Καθετήρας καθοδήγησης/Ελάχιστη εσωτερική διάμετρος
HU. Vezető katéter/Minimális belső átmérő

EN. Sterilized using ethylene oxide
DE. Sterilisiert mit Ethylenoxid
FR. Stérilisation à l'oxyde d'éthylène
IT. Sterilizzato con ossido di etilene
ET. Sterilizovano korišćenjem etilen oksida
PL. Sterylizowane przy użyciu tlenku etylenu
DA. Steriliseret ved hjælp af ethylenoxid
FI. Steriloitu etyleenioksidilla
NL. Gesteriliseerd met ethyleenoxide
CS. Sterilizováno pomocí ethylenoxidu
LV. Sterilizēts, izmantojot etilēnoksīdu
LT. Sterilizuota naudojant etileno oksidą
PT. Esterilizado com óxido de etileno
SV. Steriliseras med hjälp av etylenoxid
SK. Sterilizovaný pomocou etylénoxidu
TR. Etilen oksit kullanılarak sterilize edilmiştir
ES. Esterilizado con óxido de etileno
EL. Αποστειρώνεται με χρήση οξειδίου του αιθyleνίου
HU. Etílen-oxiddal sterilizálva

EN. Keep Dry
DE. Trocken aufbewahren
FR. Conservez-le au sec
IT. Conservare in luogo asciutto
ET. Säilitada kuivas kohas
PL. Utrzymuj w suchości
DA. Opbevarer på et tørt sted
FI. Säilytä kuivassa paikassa
NL. Droog bewaren
CS. Skladujte na suchém místě
LV. Uzglabāt sausu
LT. Laikykite sausoje vietoje
PT. Manter seco
SV. Förvaras på en torr plats
SK. Udržujte suché
TR. Kuru tut
ES. Mantener seco
EL. Αποθηκεύστε σε ξηρό μέρος
HU. Száraz helyen tárolni

EN. Catalogue Number
DE. Katalognummer
FR. Numéro du catalogue
IT. Numero di catalogo
ET. Kataloogi number
PL. Numer katalogu
DA. Katalognummer
FI. Luettelon numero
NL. Catalogusnummer
CS. Katalogové číslo
LV. Kataloga numurs
LT. Katalogo numeris
PT. Número de catálogo
SV. Katalognummer
SK. Katalógové číslo
TR. Katalog numarası
ES. Número de catálogo
EL. Αριθμός καταλόγου
HU. Katalógusszám

EN. Batch Code
DE. Chargennummer
FR. Code de lot
IT. Codice lotto
ET. Partii kood
PL. Kod partii
DA. Parti kode
FI. Erän koodi
NL. Lot code
CS. Kód pozemku
LV. Partijas kods
LT. Sklypo kodas
PT. Código de lote
SV. Kod för partiet
SK. Kód šarže
TR. Lot kodu
ES. Código de lote
EL. Κωδικός παρτίδας
HU. Tétel Kód

EN. Do Not Use If Package Is Damaged
DE. Nicht verwenden, wenn die Verpackung beschädigt ist
FR. Ne l'utilisez pas si l'emballage est endommagé
IT. Non utilizzare se la confezione è danneggiata
ET. Mitte kasutada, kui pakend on kahjustatud
PL. Nie używać, jeśli opakowanie jest uszkodzone
DA. Må ikke anvendes, hvis emballagen er beskadiget
FI. Älä käytä jos pakkaus on vahingoittunut
NL. Niet gebruiken indien verpakking beschadigd is
CS. Nepoužívejte, pokud je obal poškozený
LV. Nelietot, ja iepakojums ir bojāts
LT. Nenaudokite, jei pakuotė pažeista
PT. Não use se o pacote estiver danificado
SV. Använd inte om förpackningen är skadad.
SK. Nepoužívajte, ak je obal poškodený
TR. Ambalaj hasarlıysa kullanmayın
ES. No lo utilice si el envase está dañado
EL. Μην χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά
HU. Ne használdja, ha a csomagolás sérült

Manufacturer: Beijing Demax Medical Technology Co., Ltd.
Address: 13-7A Jingshengnansi Street, Tongzhou District, Beijing,
101102 P.R. China
Post Code: 101102
Tel.: 0086-10-59771799
Fax: 0086-10-59771883
Email: marketing@demaxmedical.com

After-sales Service: Beijing Demax Medical Technology Co., Ltd.
Address: 13-7A Jingshengnansi Street, Tongzhou District, Beijing, 101102 P.R.China
Post Code: 101102
Tel.: 0086-10-59771799
Fax: 0086-10-59771883
Email: marketing@demaxmedical.com

EU Authorized Representative:
Shanghai International Holding Corp. GmbH (Europe)
Eiffelstrasse 80, 20537 Hamburg, Germany
Tel.: 0049-40-2513175
Fax: 0049-40-255726

EN. Date Of Manufacture
DE. Datum der Herstellung
FR. Date de fabrication
IT. Data di produzione
ET. Tootmise kuup äev
PL. Data produkcji
DA. Produktionsdato
FI. Tuotantopäivä
NL. Productiedatum
CS. Datum výroby
LV. Ražošanas datums
LT. Pagaminimo data
PT. Data de fabricação
SV. Produktionsdatum
SK. Datum výroby
TR. Üretim tarihi
ES. Fecha de fabricación
EL. Ημερομηνία παραγωγής
HU. Gyártás dátuma

EN. Contents
DE. Inhalt
FR. Contenus
IT. Contenuti
ET. Sisu
PL. Zawartość
DA. Indhold
FI. Sisältö
NL. Inhoud
CS. Obsah
LV. Saturs
LT. Turinys
PT. Conteúdo
SV. Innehåll
SK. Obsah
TR. İçindekiler
ES. Contenido
EL. Περιεχόμενα
HU. Tartalomjegyzék

EN. Maximum Guidewire Diameter
DE. Maximaler Durchmesser des Führungsdrahtes
FR. Diamètre maximal du fil-guide
IT. Diametro filo guida massimo
ET. Maksimaalne juhttraadi läbimõõt
PL. Maksymalna średnica przewodu ołowianego
DA. Maksimal diameter af styret led
FI. Ohjauslangan enimmääräinen läpimittaus
NL. Maximale diameter geleidingsdraad
CS. Aximální průměr vodičového drátu
LV. Maksimālais vadotnes diametrs
LT. Maksimalus kreipiamasis laidas skersmuo
PT. Diâmetro máximo do fio-guia
SV. Maximal diameter för styrtråden
SK. Maximální priemer vodičového drátu
TR. Maksimum kurşun tel çapı
ES. Diámetro máximo del cable guía
EL. Μέγιστη διάμετρος οδηγού σύρματος
HU. Maximális vezetékvezető átmérő

EN. Medical device
DE. Medizinisches Gerät
FR. Dispositif médical
IT. Dispositivo medico
ET. Meditsiiniline seade
PL. Narzędzia medyczne
DA. Medicinsk udstyr
FI. Lääkinnälliset välineet
NL. Medische apparatuur
CS. Zdravotnické zařízení
LV. Medicīniska iekārta
LT. Medicininis prietaisas
PT. Aparelho médico
SV. Medicinsk utrustning
SK. Zdravotnícké zařízení
TR. Tıbbi cihaz
ES. Dispositivo médico
EL. Ιατρική Συσκευή
HU. Orvosi eszköz

Cobalt CAS No.
7440-48-4

EN. Contains hazardous substances
DE. Enthält gefährliche Stoffe
FR. Contient des substances dangereuses
IT. Contiene sostanze pericolose
ET. Sisaldab ohtlikke aineid
PL. Zawiera substancje niebezpieczne
DA. Indeholder farlige stoffer
FI. Sisältää vaarallisia aineita
NL. Bevat gevaarlijke stoffen
CS. Obsahuje nebezpečné látky
LV. Satur bīstamas vielas
LT. Sudėtyje yra pavojingų medžiagų
PT. Contém substâncias perigosas
SV. Innehåller farliga ämnen
SK. Obsahuje nebezpečné látky
TR. Tehlikeli maddeler içerir
ES. Contiene sustancias peligrosas
EL. Περιέχει επικίνδυνες ουσίες
HU. Veszélyes anyagokat tartalmaz

EN. Do Not Resterilize
DE. Nicht resterilisieren
FR. Ne stérilisez pas
IT. Non ristilizzare
ET. Ärge steriliseerige uuesti
PL. Nie sterylizuj ponownie
DA. Må ikke resteriliseres
FI. Ei saa steriloida
NL. Niet hersteriliseren
CS. Nesterilizujte znovu
LV. Nesterilizēt atkārtoti
LT. Nesterilizuoti pakartotinai
PT. Não reesterilize
SV. Omsteriliseras inte
SK. Neresterilizujte
TR. Yeniden sterilize etmeyin
ES. No reesterilizar
EL. Μην επαναστεριώνετε
HU. Ne sterilizálja újra

EN. Manufacturer
DE. Hersteller
FR. Fabricant
IT. Produttore
ET. Tootja
PL. Producent
DA. Producent
FI. Valmistaja
NL. Producent
CS. Výrobce
LV. Ražotājs
LT. Gamintojas
PT. Fabricante
SV. Tillverkare
SK. Výrobca
TR. Üretici
ES. Fabricante
EL. Κατασκευαστής
HU. Gyártó

EN. CE Marking Of Conformity
DE. CE Konformitätskennzeichnung
FR. Marquage CE de conformité
IT. Marchio CE di Conformità Europea
ET. CE-vastavusmärgis
PL. Znak zgodności CE
DA. CE-mærkning af overensstemmelse
FI. CE-vaatimustenmukaisuusmerkintä
NL. CE markering van conformiteit
CS. Označení shody CE
LV. CE atbilstības marķējums
LT. CE atitikties ženklas
PT. Marcação CE de conformidade
SV. CE-märkning av överensstämmelse
SK. Označenie zhody CE
TR. CE uygunluk işareti
ES. Marcado CE de conformidad
EL. Σήμανση συμμόρφωσης CE
HU. CE megfeleléségi jelölés

EN. Consult Instruction For Use
DE. Gebrauchsanweisung beachten
FR. Consultez le mode d'emploi
IT. Consultare le istruzioni per l'uso
ET. Vi kasutusjuhend
PL. Zapoznaj się z instrukcją obsługi
DA. Se brugervejledningen
FI. Katso käyttöohje
NL. Raadpleeg de gebruiksaanwijzing
CS. Viz uživatelská příručka
LV. Skatiet lietošanas instrukciju
LT. Žr. naudotojų vadovą
PT. Consulte as instruções de uso
SV. Se bruksanvisningen
SK. Pozrite si návod na obsluhu
TR. Çalıştırma talimatlarına bakın
ES. Consultar las instrucciones de uso
EL. Ανατρέξτε στο εγχειρίδιο χρήστη
HU. Lásd a felhasználói kézikönyvet

EN. Unique device identifier
DE. Eindeutige Geräteerkennung
FR. Identifiant unique du dispositif
IT. Identificatore unico del dispositivo
ET. Unikaalne seadme identifikaator
PL. Unikalny identyfikator urządzenia
DA. Unik enhedsidentifikation
FI. Yksilöllinen laitetunniste
NL. Unieke apparaatidentificatie
CS. Jedinečný identifikátor zařízení
LV. Unikāls ierīces identifikators
LT. Unikalus įrenginio identifikatorius
PT. Identificador exclusivo do dispositivo
SV. Unik enhetskod
SK. Jedinečný identifikátor zariadenia
TR. Benzersiz cihaz tanımlayıcısı
ES. Identificador único del dispositivo
EL. Μοναδικό αναγνωριστικό συσκευής
HU. Egyedi eszközazonosító

EN. Nominal Pressure
DE. Berstdruck
FR. Pression d'éclatement
IT. Pressione nominale di scoppio
ET. Lõhkemisrõhk
PL. Ciśnienie rozrywające
DA. Tryk ved sprængning
FI. Murtumispaine
NL. Barstdruk
CS. Tlak při výbuchu
LV. Sprādziena spiediens
LT. Sprogimo slėgis
PT. Pressão de ruptura
SV. Tryck vid utbrott
SK. Trhac tlak
TR. Patlama basıncı
ES. Presión de rotura
EL. Πίεση έκρηξης
HU. Berobbanási nyomás

EN. Nominal Pressure
DE. Nenndruck
FR. Pression nominale
IT. Pressione nominale
ET. Nominaalne rõhk
PL. Ciśnienie nominalne
DA. Nominelt tryk
FI. Nimellinen paine
NL. Nominale druk
CS. Jmenovitý tlak
LV. Nominālais spiediens
LT. Nominalusis slėgis
PT. Pressão nominal
SV. Nominellt tryck
SK. Nominální tlak
TR. Nominal basınç
ES. Presión nominal
EL. Ονομαστική πίεση
HU. Névleges nyomás

EN. Keep Away From Sunlight
DE. Von Sonnenlicht fernhalten
FR. Tenez-le à l'abri des rayons de soleil
IT. Tenere lontano dalla luce del sole
ET. Hoida eemal päikesevalgusest
PL. Trzymaj się z dala od światła słonecznego
DA. Holdes væk fra sollys
FI. Pidäpoissa auringonvalolta
NL. Blijf uit de buurt van zonlicht
CS. Chraňte před slunečním světlem
LV. Paliecties prom no saules gaismas
LT. Saugokite nuo saulės spindulių
PT. Mantenha longe da luz do sol
SV. Håll dig borta från solljus
SK. Drž sa dale od slnečného svetla
TR. Güneş ışığından uzak durun
ES. Mantener alejado de la luz solar
EL. Μακριά από το ηλιακό φως
HU. Tartsa távol a napfénytől

EN. Single Sterile Barrier System With Protective Packaging Inside & Sterilized Using Ethylene Oxide
DE. Einzelnes steriles Barriersystem mit Schutzverpackung im Inneren & sterilisiert mit Ethylenoxid
FR. Système de barrière stérile unique avec emballage de protection à l'intérieur et stérilisé à l'oxyde d'éthylène
IT. Sistema di barriera sterile doppia all'interno della confezione protettiva e sterilizzata usando ossido di etilene
ET. Sterilne tihokordse barjääri süsteem koos sisepakendiga ja steriliseeritud etüleenoksiidiga
PL. Sterylny system z pojedynczą barierą z wewnętrznym opakowaniem i sterylizowany tlenkiem etylenu
DA. Sterilt enkeltbarriersystem med indvendig emballage og steriliseret med ethylenoxid
FI. Steriili yksisuulkujärjestelmä jossa on sisäpakkaus ja joka on steriloitu etyleenioksidilla.
NL. Steriel enkelvoudig barrièrestelsysteem met binnenverpakking en gesteriliseerd met ethyleenoxide
CS. Sterilní jednobariérový systém s vnitřním obalem a sterilizovaný ethylenoxidem
LV. Sterila viena barjera sistēma ar iekšējo iepakojumu un sterilizēta ar etilēna oksīdu
LT. Sterili vieno barjero sistema su vidine pakuote, sterilizuota etileno oksidu
PT. Sistema de barreira estéril única com embalagem protetora interna e esterilizada usando óxido de etileno
SV. Sterilt enkelbarriärsystem med innerförpackning som steriliseras med etylenoxid.
SK. Sterilní jednobariérový systém s vnitřním obalem a sterilizovaný etylénoxidom
TR. İç ambalajlı ve etilen oksit ile sterilize edilmiş steril tek bariyer sistemi
ES. Sistema de barrera estéril único con embalaje protector interior y esterilizado utilizando óxido de etileno
EL. Αποστειρωμένο σύστημα ενιαίου φραγμού με εσωτερική συσκευασία και αποστείρωση με οξείδιο του αιθυλενίου
HU. Steril egysoros barrier rendszer belső csomagolással és etilén-oxidallal sterilizálva

ENGLISH

CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE. OBSERVE ALL WARNINGS AND PRECAUTIONS NOTED THROUGHOUT THESE INSTRUCTIONS. FAILURE TO DO SO MAY RESULT IN COMPLICATIONS.

Description

The PTCA Balloon Dilatation Catheter (Gusta® /Gusta® NC/ Gusta® II/ Gusta® II NC) is designed to allow rapid exchange of the catheter using a standard length guidewire. The proximal shaft is a PTFE coated stainless steel tube. The steel construction allows superior proximal pushability with a smooth transition to a distal shaft specifically designed to be traceable through tortuous vasculature. The distal portion of the catheter is coated with a hydrophilic coating. It is composed of polyvinyl pyrrolidone and polyurethane. No contraindications related to the coating was found. The balloon material allows high-pressure dilatation while maintaining control of balloon size. Two radiopaque platinum marker bands are located within the balloon with the exception of 1.25 mm and 1.5 mm diameter balloons, where the marker band is centrally located. The catheter is compatible with a 1.47 mm(5F) or larger guiding catheter. The internal lumen of the catheter accepts a standard 0.36 mm (0.014 inch) PTCA guidewire. The proximal part of the guidewire enters the catheter's tip, and advances coaxially out the catheter's proximal port, thereby allowing both coaxial guidance and rapid exchange of catheters with a single standard length guidewire. Two marked sections are located on the proximal shaft to indicate catheter position relative to the tip of either a brachial or femoral guiding catheter. The design of this dilatation catheter does not incorporate a lumen for distal dye injections and distal pressure measurements.

Trade-mark	series	Compliance
Gusta®	DBC	semi-compliant
Gusta®NC	DNC	non-compliant
Gusta®II	DPC	semi-compliant
Gusta®IINC	DHC	non-compliant

Information of Basic UDI-DI:

The SSCP is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI.

URL to the Eudamed public website: <https://ec.europa.eu/tools/eudamed>
Basic UDI-DI:69437182GustaDX

Intended purpose

The PTCA Balloon Dilatation Catheter is intended for dilatation of the stenotic portion (non-CTO) of coronary arteries for the purpose of improving myocardial perfusion.

Intended Users

The PTCA Balloon Dilatation Catheter is meant to be used by experienced physicians trained in coronary angiography and/or percutaneous transluminal coronary angioplasty (PTCA) procedures. It is solely the hospital/clinic's responsibility to assign qualified users to use the PTCA Balloon Dilatation Catheters.

Patient Target Group

The DBC series DPC series PTCA semi-compliant Balloon Dilatation Catheter are intended for patients suffering from non-CTO coronary artery disease (CAD) requiring percutaneous transluminal coronary angioplasty (PTCA) procedure (pre-dilatation) before further percutaneous coronary interventions (PCI).

The DNC series) and DHC series PTCA compliant Balloon Dilatation Catheters are intended for patients suffering from non-CTO coronary artery disease (CAD) requiring either percutaneous transluminal coronary angioplasty (PTCA) procedure (pre-dilatation) before further percutaneous coronary interventions (PCI) or requiring in-stent dilatation or both.

How Supplied

Contents:

- One (1) PTCA Balloon Dilatation Catheter
- One (1) Flush needle

Sterile: Sterilized with ethylene oxide gas. Non-pyrogenic.

Storage: Store in a dry, dark, cool place.

Expired Date: 3 years after ETO sterilization.

Transportation condition:

keep away from heavy and sunshine, keep dry

Temperature limit: -18°C~55°C

Humidity limitation: 0~85%

Indication for use

The PTCA Balloon Dilatation Catheter (Gusta® DBC series/Gusta® II DPC series) is indicated

for pre-dilatation of the stenotic portion (non-CTO) of a coronary artery during percutaneous transluminal coronary angioplasty (PTCA)

procedure.

The PTCA Balloon Dilatation Catheter (Gusta® NC DNC series/Gusta® II NC DHC series) is indicated for pre-dilatation of the stenotic portion (non-CTO) of a coronary artery during percutaneous transluminal coronary angioplasty (PTCA) procedure; and it's also indicated for post- intra stent dilatation of the stenotic portion (non-CTO) of coronary artery during percutaneous coronary intervention.

Clinical Benefit

In the recent years, percutaneous coronary interventions (PCI) have become more frequently used in the treatment of CAD, when compared with medical treatment or CABG. Pre- and post-dilatation of stenotic lesions were shown to significantly relieve anginal symptoms, resulting from an improved blood flow and dynamics.

Clinical benefits also include an improvement of the overall myocardial perfusion (TIMI 3 achieved post-procedure), decrease of thrombosis and MACE rates.

The device under evaluation can achieve such clinical benefits by reducing the stenotic portion of the coronary artery and improve/restore the overall myocardial perfusion.

Contraindications

- Unprotected left main coronary artery.
- Coronary artery spasm in the absence of a significant stenosis.

Warnings

• This device is intended for single use only. DO NOT re-sterilize and/or reuse the device, as this can potentially increase the risk of contamination by bloodborne pathogens causing infection or cross-infection. Multiple sterilization cycles or reuse may decrease the device performance.

• PTCA in patients who are not acceptable candidates for coronary artery bypass graft surgery requires careful consideration, including possible hemodynamic support during PTCA, as treatment of this patient population carries special risk.

• Use only the recommended balloon inflation medium, such as contrast agent. Never use air or any gaseous medium to inflate the balloon.

• When the catheter is exposed to the vascular system, it should be manipulated while under high quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding. Continuing to advance or retract the catheter while under resistance may result in damage to the vessels and / or damage / separation of the catheter.

• Do not use, or attempt to straighten a

catheter if the shaft has become bent or kinked; this may break the shaft. Instead, prepare a new catheter.

- Balloon pressure should not exceed the rated burst pressure (RBP) indicated on the package. The RBP is based on the results of in-vitro testing. At least 99.9% of the balloons (with a 95% confidence) will not burst at or below their RBP.
- The catheter should only be used with an appropriate device for balloon inflation that includes a pressure gauge or display. If the inflation device fails, it should be replaced.
- To reduce the potential for vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.
- PTCA should only be performed at hospitals where emergency coronary artery bypass graft surgery can be quickly performed in the event of a potentially injurious or life-threatening complication.
- Treatment of moderately or heavily calcified lesions is considered moderate risk, increasing the risk for acute closure, vessel trauma, balloon burst, balloon entrapment and associated complications. If resistance is felt, determine the cause before proceeding. Continuing to advance or retract the catheter while under resistance may result in damage to the vessels and / or damage / separation of the catheter.
- Care should be taken in handling of the distal end of the catheter as damage to the coating may cause the release of particles resulting in embolism.
- After retraction of the balloon catheter from the guiding catheter, the physician should confirm the integrity of the balloon catheter.
- In the event of catheter damage or component separation, the physician should remove the separated catheter components from the patient according to their individual condition and appropriate treatment protocol(s).
- Use the catheter prior to the "Use By" date

specified on the package.

- Inspect all products prior to use. Do not use if the package has been opened or damaged and prepare a new catheter.
- Potential biological hazards exist after use. Please dispose of product according to medical practice, local regulations and laws.
- To prevent the hazard of the sharp Flush Needle, please insert the flush needle into the flush needle cap at disposal of the device.

Precautions

- Before angioplasty, the dilatation catheter should be examined to verify functionality and ensure its size is suitable for the specific procedure for which it is to be used.
- The catheter system should be used only by physicians trained in the performance of percutaneous transluminal coronary angioplasty (PTCA).
- During the procedure, appropriate anti-coagulant and coronary vasodilator therapy must be provided to the patient as needed. After the procedure, anti-coagulant therapy should be continued for a period of time as determined by the physician.
- Percutaneous Coronary balloon angioplasty may cause temporary patient discomfort.
- CMR(Cobalt>0.10%w/w) substances incorporated into the device could result in sensitization or an allergic by the patient or user.

Adverse Event

- Cardiac death
- Target lesion revascularization
- Non-fatal myocardial infarction or acute myocardial infarction
- Stroke
- MACE event (Cardiac death, Target lesion revascularization, Non-fatal myocardial infarction, stroke)
- Intra-coronary Embolization (Air embolism)
- Coronary artery wall dissection

- Restenosis of the dilated vessel
- Intra-coronary thrombosis, acute thrombosis
- Coronary artery perforation or rupture
- Haemorrhage or hematoma
- Total occlusion of the coronary artery, artery closure
- Pseudoaneurysm.
- Arteriovenous fistula
- Arrhythmias, including ventricular fibrillation, cardiogenic shock
- Drug reaction-allergic reaction to contrast medium
- Hypo- / hypertension
- Infection
- Coronary artery spasm
- Urgent CABG surgery (emergency bypass graft surgery)
- Unstable angina

Note:

any serious incident that has occurred in relation to the device or its procedure should be reported to the manufacturer and the competent authority of the member state in which the user and/or patient is established.

Materials To Be Used In Combination With A Balloon Catheter Include:

- Suitable Guide wire, refer to label claim.
- Suitable guiding catheter, refer to label claim.
- Appropriate inflation medium.
- Pressure-indicating inflation device.
- Syringe

Performance Characteristics

- The effective length of the catheter (mm): DBC/DNC 1410 mm DHC/DPC 1390 mm
- The position of marker bands: the double marker bands are flush with the balloon shoulders a single marker band is in the center of the effective length of the balloon

- The diameter of the distal tip where the guidewire is inserted:0.016 inch
- Nominal Pressure (NP): DBC and DPC: 6 atm DNC and DHC: 12 atm
- Rated Burst Pressure (RBP): DBC and DPC: 16 atm (expect 1.25 and 1.5 diameter is14atm) DNC and DHC: 20 atm
- Maximum outer diameter of guidewire used: 0.356 mm(0.014") Minimum inner diameter of guiding catheter: 1.47 mm(5F)
- Hydrophilic Coating: The catheter shaft and balloon are coated with a hydrophilic coating to help reduce the frictional force between the catheter and the guiding catheter as well as the vessel wall.

Preparation For Use

Prior to use examine all equipment carefully for defects. Examine the dilatation catheter for bends, kinks or other damage. Do not use any defective equipment. Prepare equipment to be used following manufacturer's instructions or standard procedure. Complete the following steps to prepare the PTCA Balloon Dilatation Catheter for use:

1. Remove the guidewire stylet from the catheter tip.
2. Remove the protection sheath off the balloon.
3. Flush the PTCA Balloon Dilatation Catheter and guidewire port by standard saline solution with the flush needle.
4. Gently insert the tip of a syringe filled with heparinized normal saline over the tip of the catheter and gently inject heparinized saline into the lumen until fluid is seen exiting the guidewire Rx port.
5. Prepare an inflation device with the recommended contrast medium according to the manufacturer's instructions.
6. Evacuate air from the balloon segment

using the following procedure:

- ① Fill a 20 cc syringe or the inflation device with approximately 4 cc of the recommended contrast medium.
- ② After attaching the syringe or inflation device to the balloon inflation lumen, orient the dilatation catheter with the distal tip and the balloon pointing in a downward vertical position.
- ③ Apply negative pressure and aspirate for 15 seconds. Slowly release the pressure to neutral, allowing contrast to fill the shaft of the dilatation catheter.
- ④ Disconnect the syringe or inflation device from the inflation port of the dilatation catheter.
- ⑤ Remove all air from the syringe or inflation device barrel. Reconnect the syringe or inflation device to the inflation port of the dilatation catheter. Maintain negative pressure on the balloon until air no longer returns to the device.
- ⑥ Slowly release the device pressure to neutral.
- ⑦ Disconnect the 20 cc syringe (if used) and connect the inflation device to the inflation port of the dilatation catheter without introducing air into the system.

CAUTION: All air must be removed from the balloon and displaced with contrast medium prior to inserting the dilatation catheter into the body. Otherwise complications may occur.

Instructions For Use

1. Insert a guidewire through the hemostatic valve following the manufacturer's instructions.
2. Advance the guidewire carefully into and through the guiding catheter. Withdraw the guidewire introducer, if used.
3. Attach a torque device to the guidewire, if desired. Under fluoroscopy, proceed with accepted PTCA techniques to advance the guidewire to and across the lesion.
4. Backload the distal tip of the dilatation catheter onto the guidewire ensuring that the guidewire exits the catheter at approximately 25 cm proximal to the balloon.
5. Advance the dilatation catheter over the

guidewire until it approaches the hemostatic valve.

6. Open the hemostatic valve. Insert the dilatation catheter while maintaining the guidewire position and tighten the hemostatic valve. To facilitate insertion, the balloon must be fully deflated to negative pressure.
7. Tighten the hemostatic valve to create a seal around the dilatation catheter without inhibiting movement of the dilatation catheter. This will allow continuous recording of proximal coronary artery pressure.
NOTE: It is important that the hemostatic valve be closed tightly enough to prevent blood leakage around the dilatation catheter shaft, yet not so tight that it restricts the flow of contrast into and out of the balloon or restricts guidewire movement.

8. Advance the dilatation catheter until the appropriate proximal marker aligns with the hemostatic valve hub. This indicates that the dilatation catheter tip has reached the guiding catheter tip.

9. Advance the dilatation catheter over the guidewire and into the stenosis. Continue under fluoroscopy and use the radiopaque marker band(s) to position the balloon's usable (dilating) section within the stenosis.

NOTE: When using the dual wire technique, a dual hemostatic valve should be used and care should be taken when introducing, torquing and removing one or both wires to avoid entanglement. During the dual wire procedure, guidewires should not be rotated more than 180 degrees in either direction, during the dual wire procedure. It is recommended that one wire be completely withdrawn from the patient before removing additional equipment.

10. Inflate the balloon to dilate the lesion using standard PTCA techniques.

11. After each subsequent inflation, the distal blood flow should be assessed.

12. If a significant stenosis persists, successive inflations may be required to resolve the stenosis. Do not exceed the rated burst pressure.

13. Apply negative pressure to completely

deflate the balloon. Fluoroscopically verify that the balloon has been deflated.

14. Confirm the results with fluoroscopy.

Exchange Procedure Technique

The PTCA Balloon Dilatation Catheter has been designed for rapid, single operator balloon exchanges. To perform a dilatation catheter exchange:

1. Loosen the hemostatic valve.
2. Hold the guidewire and hemostatic valve in one hand, while grasping the balloon shaft in the other hand.
3. Maintain the guidewire position in the coronary artery by holding the wire stationary, and pulling the dilatation catheter out of the guiding catheter while monitoring the wire position under fluoroscopy.
4. Withdraw the deflated dilatation catheter until the guidewire lumen is reached. Carefully pull back the flexible, distal portion of the dilatation catheter out of the rotating hemostatic valve while maintaining the guidewire's position across the lesion.
5. Slide the distal tip of the dilatation catheter out of the hemostatic valve, and tighten the valve onto the guidewire to hold it securely in place.
6. Prepare the next dilatation catheter to be used, as previously described in the Preparation for the Use section.
7. Reload another dilatation catheter onto the guidewire as previously described under the Instructions for Use Section, Step 4, and continue the procedure accordingly.

eIFU

The exact same pdf version e-IFU can also be found on the website of Demax Medical: https://www.demax.group/news_detail/57.html.

You can readable using adobe or other PDF reader program.

If you cannot download it from the website, please contact the manufacturer:

Tel.: 0086-10-59771799

Fax: 0086-10-59771883

Note: When the manufacturer's instruction for use is updated, it will be uploaded timely. For it is difficult to trace to every end user to inform the change, so we advise the customer to browse and check it regularly.

References

The physician should consult recent literature on current medical practice on balloon dilatation, such as that published by ACC/AHA or ESC/EAS.

Disclaimer Of Warranty And Limitation Of Remedy

Descriptions or specifications in Demax Medical printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. Demax Medical will not be responsible for any direct, incidental, or consequential damages resulting from the misuse of the product.

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